


FS-200L

Single-Use Surgical Scissors

FS-410U

Single-Use Loop Cutter



Single-Use Surgical Scissors FS-200



Single-Use Loop Cutter FS-410

Introducing single-use surgical scissors and loop cutters

Product Features

■ FS-200L

The scissors are precisely designed and manufactured to realize smooth incision of tissue without high frequency power supply.





■ FS-410U

With their characteristic one side open, the scissors are precisely designed and manufactured to realize smooth cutting of the residual end of the loop used for ligating tissue.



Specifications

Article Number	Model	Quantity	Working Length	Min. Channel Size
N5781230	FS-200L	1/box	1,619 mm	2.8 mm
N5781130	FS-410U	1/box	2,300 mm	2.8 mm



Olympus reserves the right of errors, modification and changes of the service and/or product offerings.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com