


X-TRAC

Extraction Bag for Foreign Bodies and Polyps



Features and Benefits

• Strong Solid Loop

Facilitates catching the object without the sides of the bag sticking together

• No Sticking of Mucosa in the Bag

Very easy transfer into histology box

• Flexible Bag

Bag can collect small, big, and even multiple pieces at a time

• Protection from Sharp Objects


Optimal isolation of sharp foreign bodies

• Specially Perforated Bag

Perforation allows tissue liquid to pass out of the bag

Please refer to the Olympus ET catalogue for more information

www.olympus.eu/et-catalogue



X-TRAC Specifications and Indications

Article Name	Article Code	Quantity/ Pack	Bag Ø mm	Min. Working Channel Ø mm	Working Length (mm)	Type of Bag	Indication	Object to Remove
Extraction Bag X-TRAC, 25 mm	E0420200	5	25	2.8	2300	Non-perforated	Gastroenterology Colonoscopy Bronchoscopy	Foreign bodies
Extraction Bag X-TRAC, 35 mm	E0420201	5	35	2.8	2300	Non-perforated	Gastroenterology Colonoscopy Bronchoscopy	Foreign bodies
Extraction Bag X-TRAC, perforated, 25 mm	E0420202	5	25	2.8	2300	Perforated	Gastroenterology Colonoscopy	Polyps
Extraction Bag X-TRAC, perforated, 35 mm	E0420203	5	35	2.8	2300	Perforated	Gastroenterology Colonoscopy	Polyps
Extraction Bag X-TRAC, 15 mm	E0420204	5	15	2.0	1600	Non-perforated	Bronchoscopy	Foreign bodies, blood, and mucus
Extraction Bag X-TRAC, 10 mm	E0420205	5	10	1.2	1200	Non-perforated	Bronchoscopy	Foreign bodies, blood, and mucus

Manufactured by:

Carmonja® GmbH
Am Wegmannbühl 4 | 83246 Unterwössen | Germany

Distributed by:

Olympus Europa SE & Co. KG
Wendenstrasse 14-18 | 20097 Hamburg | Germany

Olympus reserves the right of errors, modification and changes of the service and/or product offerings.