

OLYMPUS[®]

Your Vision, Our Future

Endo**Therapy**

Single Use Spray Catheter PW-205V

Efficient Spraying During Endoscopy



Maximize the Effectiveness and Efficiency of Spraying with the Single Use Spray Catheter PW-205V from Olympus

Main Features

■ Wider, more evenly distributed spraying

With the PW-205V, you can spray dye, saline, or another solution more evenly and over a wider area than with a syringe attached to the endoscope's biopsy port. These features can help shorten procedure time and ease patient discomfort.

■ Kink-resistant design

The stylet provided with the PW-205V ensures smooth insertion into the endoscope and prevent the catheter from kinking.

■ Single-use design

The PW-205V's single-use design optimizes device handling for each endoscopic procedure.



Specifications

Model	PW-205V
Shape of the Distal End	Spray type (with nozzle)
Maximum Insertion Portion Diameter	2.45 mm
Working Length	2400 mm
Minimum Channel Size	2.8 mm



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of
sales and distribution locations visit:
www.olympus.com